


# C, S determination in Ferrochrome (Low Carbon)



**Suitable analyzers**

- ELEMENTRAC CS-*i*

**Used accessories**

- Ceramic crucibles (90149)
- Tungsten (90220)
- High purity iron accelerator (88600-0013)
- Suitable calibration material (NIST or other)



**Application Settings**

**I) General**

Sample type: Advanced  
 Standby flow: 10 l/h  
 Purging while closing: yes  
 Open Furnace: yes  
 Furnace purge through: Exhaust  
 Furnace purge time: 3 sec  
 Furnace purge flow: 180 l/h

**Stabilizing**

Lance valve: on  
 Stabilize by time: off  
 Stability: 0.02 V  
 Minimum time: 15 sec  
 Maximum time: 30 sec

**II) Analysis**

Voltage: 100 %  
 Power duration: 35 sec  
 Flow: 180 l/h  
 Chamber only: 1 sec  
 Lance and chamber: 1 sec  
 Drift compensation: on


Channel	Max time [sec]	Min time [sec]	Integration delay [sec]	Comparator factor [%]
Low C	80	30	7	0.1
Low S	80	15	10	0.1

**III) Postwaiting**

Postwaiting time: 10 sec

Subject to technical modification and errors

## C, S determination in Ferrochrome (Low Carbon)



### Sample preparation

Make sure that your sample is free from contaminations and inclusions which could influence the carbon determination. Pre-heat the crucibles at least for 1 h at 1000 °C. Let the crucibles cool down in a desiccator.

### Procedure

- Prepare ELTRA analyzer (e.g. exchange anhydron, sodium hydroxide, platin catalyst if necessary); clean the combustion tube, brush, heat shield, dust trap
- Run three warm up samples (e.g. steel samples (92400-3050) with a minimum weight of 500 mg; add 1.7 g tungsten)
- Calibrate the analyzer with suitable calibration material (NIST or other)

The procedure of analysis Ferro-Chromium should be like this:

- (1) Weigh in approx. 200 mg of sample into the crucible
- (2) Add 0.7 g of high purity iron accelerator (88600-0013)
- (3) Add 1.7 g of tungsten (90220)


Repeat steps (1) – (3) at least three times;

Mark the results and use the calibration function in the software.

**-> Now start with the actual analysis.**



Typical results		
Euronorm 580-1 <sup>1)</sup>		
Weight (mg)	Carbon (ppm)	Sulfur (ppm)
198.4	190.62	22.75
202.6	181.32	19.21
208.3	186.46	17.45
202.4	181.99	17.68
198.3	189.73	16.45
200.7	197.63	16.86
200.8	186.47	18.37
201.5	198.48	16.66
204.9	199.25	16.40
199.6	190.05	16.62
Mean value		
	190.20	17.84
Deviation / Relative deviation (%)		
	6.50/3.42	1.95/10.94




<sup>1)</sup> certified value: C: 0.019% ±0.002 (10.5%)  
S: not certified

# C, S determination in Ferrochrome (Low Carbon)



Typical results		
D-Lab (Sweden): DFS 1 <sup>1)</sup>		
Weight (mg)	Carbon (ppm)	Sulfur (ppm)
200.5	563.23	17.26
199.2	547.90	19.95
203.3	551.84	16.67
200.5	549.86	18.28
200.1	544.32	18.05
204.6	540.70	17.62
200.3	546.44	20.31
201.4	538.57	16.54
202.8	566.90	14.05
202	541.66	13.82
Meanvalue		
	549.14	17.25
Deviation / Relative deviation (%)		
	9.39/1.71	2.14/12.42



<sup>1)</sup> certified value: C: 0.0549 ±0.0014 (2.56%)  
S: 0.0017% (Deviation not certified)