


C/S determination in limestone



Suitable analyzers

- ELEMENTRAC CS-d (Resistance Furnace)

Used accessories

- Disposable porcelain boats (90160)
- Suitable calibration material



Application Settings


I) General

Temperature:	1450 °C (±20 °C tolerance)
Stabilize baseline:	Enable
Stability:	0.01 V
Minimum time:	20 sec
Maximum time:	60 sec
Flow:	180 l/h
Drift compensation:	Disable

II) Analysis

Channel	Min time [sec]	Max time [sec]	Integration delay [sec]	Peak max [V]	Comperator level [mv]	Comperator peak [%]
High Carbon	90	300	5	8	20	1
Low Sulfur	200	300	5	8	5	1

C/S determination in limestone



Sample preparation

Dry the sample to constant mass at 105 °C (at least 1 hour).

Procedure

- Prepare and clean the ELTRA analyzer (e.g. exchange anhydron, filter, boat stop) and set the furnace temperature to 1450 °C
- Run at least three warm up samples (e.g. ELTRA 92511-3020) with a medium sample weight of 250 mg until the results are consistent
- Calibrate the system with a suitable calibration material (NIST or other):

- (1) Weigh in 250 mg of sample in a porcelain boat (90160)
- (2) Start analysis (F5 Button)
- (3) Wait until baseline is stable
(Look at message in ELEMENTS software and wait for green light at the resistance furnace)
- (4) Load the sample into the furnace and wait until the PC calculates results
- (5) Remove combustion boat


Repeat steps (1) – (5) at least three times;

Mark the results and use the calibration function in the software.

-> Now start with the actual analysis.



Typical results		
ELTRA 90812-3002 (LOT 101602) *1		
Weight (mg)	Carbon (%)	Sulfur (%)
258.6	11.6	0.44
252.6	11.6	0.43
271.6	11.6	0.43
269.5	11.9	0.41
280.6	11.6	0.42
251.7	11.6	0.40
269.7	11.6	0.40
275.7	11.7	0.40
273.2	11.7	0.40
289.7	11.7	0.40
Average Values		
	11.70	0.41
Deviation / Relative deviation (%)		
	0.09 (0.8%)	0.01 (3.4%)



*1 certified values: C: 11.72 % ±0.40 ; S: 0.418 % ±0.04