


C, S determination in ore samples



Suitable analyzers

- CS-800
- CS-2000 (induction furnace)

Used accessories

- Ceramic crucibles (90150)
- Tungsten (90220)
- Iron accelerator (90260)
- Suitable calibration material (NIST or other)

Settings

- | | |
|----------------------------|---------------------------------|
| ■ Comparator level: 20 mV | ■ Base line deviation: 10 mV |
| ■ Minimum time: 40 sec | Base line time: 15 sec |
| ■ Maximum time: 3:00 min | ■ T1 - Pre-purging: 5 sec |
| ■ Integration delay: 6 sec | ■ T2 - Chamber only: 5 sec |
| ■ Post-waiting: 15 sec | ■ T3 - Chamber and lance: 3 sec |
| ■ Pre-purging: 10 sec | ■ T4 - Generator: 45 sec |



CS-800

CS-2000
(induction furnace)

Sample preparation

For best results grind the sample down to a particle size of approx. 200 µm.
 Dry the sample to constant mass at 105°C (at least 1 hour).
 Pre-heating the crucibles (1 h; 1000°C) could improve the precision.

Procedure

- Prepare ELTRA analyzer (e.g. exchange anhydrous, sodium hydroxide, copper oxide when necessary); clean the combustions tube
- Run three warm up samples (e.g. steel samples (92400-3050) with a minimum weight of 500 mg; add 1.5 g tungsten)
- Calibrate the analyzer with suitable calibration material (NIST or other)
 - (1) Weigh in 150 mg of sample in the ceramic crucibles (90150)
 - (2) Add 0.7 (±0.1) g of iron accelerator (90260) and 2.0 (±0.1) g of tungsten (90220)
 - (3) Place the crucible on the pedestal and start analysis

Repeat Step (1) – (3) at least three times;

Mark the results and use the calibration function in the software.

-> Now start with the actual analysis.

Typical results

AR 4006 : Iron ore

Weight (mg)	% C	% S
162.8	4.18	4.14
172.3	4.13	4.12
196.4	4.14	3.94
159.8	4.16	4.14
163.1	4.21	3.97
169.2	4.20	4.04
174.3	4.24	4.10
168.4	4.24	4.04
167.6	4.23	4.14
170.9	4.17	4.08

Average values

	4.19	4.07
--	------	------

Deviation

	0.04 / 0.95 %	0.07 / 1.77 %
--	------------------	------------------