

Oxygen and nitrogen determination in steel samples



Suitable analyzers

- ELEMENTRAC ONH-*p*
- ELEMENTRAC ON-*p*

Used accessories

- Graphite crucibles (90180 and 90185)
- Suitable calibration material (NIST or other)



Application Settings

I) General

Furnace mode: ON

Furnace cooling: 35/45 °C

Standby Flow: 0

A flow of 10 l/h could improve precision when there is a long time distance between 2 measurements.

II) Outgasing and stabilizing

Setting / Phase	Time [sec]	Power [W]	Flow [l/h]
Outgasing	60	4600	27
Stabilizing	60	4400	27

A second outgasing cycle or an increased outgasing time could improve the precision for very low oxygen and nitrogen contents.

III) Analysis

Power duration: 80 sec

Drift compensation: on

Power: 4400 W

Open furnace: yes

Flow: 27 l/h

Channel	Minimum time [sec]	Maximum time [sec]	Integration delay [sec]	Comparator factor [%]
Low and High O	35	80	5	0.5
Low and High N	35	80	10	0.5

IV) Postwaiting

Postwaiting time: 15 sec

Furnace clean up: No

Oxygen and nitrogen determination in steel samples



Sample preparation

Make sure that the surface of the steel is free from contaminations; otherwise clean the sample with acetone p.a. and let dry at atmosphere.

Procedure

- Prepare ELTRA analyzer (exchange anhydron, sodium hydroxide, Schuetze reagent if necessary), clean furnace, sample drop mechanism, electrode tip (if necessary)
- Run three blanks with empty crucibles
- Calibrate the analyzer with suitable calibration material (NIST or other)
 - (1) Fill one empty inner crucible (90180) in one outer crucible (90185) and place them on the electrode tip, close furnace
 - (2) Weigh calibration material and place it in the sample drop mechanism and start analysis
 - (3) After analysis give the inner crucible into waste and fill in a new one. The outer crucibles can be used approximately 10 times

Repeat steps (1) – (3) at least three times;
 Mark the results and use the calibration function in the software.

-> Now start with the actual analysis.



Typical results		
Steel ELTRA 91100-1001 (LOT 714 A)		
Weight (mg)	ppm O	ppm N
994.9	6.1	18.8
999.0	5.1	18.0
1000.2	5.7	17.8
997.6	6.3	18.4
1000.4	6.9	18.6
997.4	5.9	17.8
997.5	7.1	19.4
994.7	5.6	18.4
996.9	5.6	19.4
998.1	5.8	18.4
Average values		
	6.0	18.5
Deviation / Relative deviation (%)		
	0.6 / 10.3 %	0.6 / 3.1 %

Typical results		
Steel ELTRA 91100-1004 (LOT 914 F)		
Weight (mg)	ppm O	ppm N
1005.4	108.3	54.8
994.3	116.2	56.1
970.4	112.2	53.7
994.6	116.4	56.6
997.5	113.6	56.1
1000.9	111.6	56.1
1000.3	114.4	54.1
1001.4	108.4	56.2
1004	107	56.8
996.1	111.8	55.1
Average values		
	111.99	55.56
Deviation / Relative deviation (%)		
	3.29 / 2.93 %	1.07 / 1.92 %