



**HOLDER FOR DOSING NOZZLES**

The holder for dosing nozzles for manual operation can be positioned anywhere on the device. For automatic operation with the Rubin 500 the nozzles are fastened directly to the grinding head.



**Topas 130**

The dosing unit for diamond suspensions and lubricants can be used universally.

The dosing interval is adjustable.

In automatic mode (Rubin 500), the dosage is automatically switched on and off in the process.

A reverse rinsing function keeps the hoses clean.

**FUNCTION WITH SAPHIR DEVICES**

SAPHIR	manual mode	automatic mode
320		
330		
350		
360		
520/Rubin 500		
530/Rubin 500		
550/Rubin 520	*	
550/Rubin 540		
560/Rubin 520	*	

\* recommended only without dosing unit for Rubín 520



BASIC MODULE

TOPAS 130



Order No.: M6300040

- » automatic dosing unit (also other brands)
- » manual or automatic function
- » variable dosing, reverse rinsing function

Diamond suspension	3x
Lubricant	1x
WxHxD	180x240x195 mm
Weight	2 kg

EQUIPMENT



POWER SUPPLY

Equipment 1  
230 V/50 Hz (1Ph/N/PE)

Order No.: A6300001

Equipment 2  
110 V/60 Hz (1Ph/N/PE)

Order No.: A6300002

HOLDER FOR DOSING NOZZLES



Equipment 1  
**automatic mode**  
Order No.: A6300003



Equipment 2  
**manual mode**  
Order No.: A6300004

